CARDIFF UNIVERSITY PRIFYSGOL CAERDYD

**ORCA – Online Research @ Cardiff** 

## This is an Open Access document downloaded from ORCA, Cardiff University's institutional repository:https://orca.cardiff.ac.uk/id/eprint/111982/

This is the author's version of a work that was submitted to / accepted for publication.

## Citation for final published version:

Baxani, Divesh K., Jamieson, William D., Barrow, David A. and Castell, Oliver K. 2018. An encapsulated droplet interface bilayer array for the high-throughput optical measurement of lipid membranes with single bilayer resolution. Biophysical Journal 114 (3), 686A. 10.1016/j.bpj.2017.11.3699

Publishers page: https://doi.org/10.1016/j.bpj.2017.11.3699

## Please note:

Changes made as a result of publishing processes such as copy-editing, formatting and page numbers may not be reflected in this version. For the definitive version of this publication, please refer to the published source. You are advised to consult the publisher's version if you wish to cite this paper.

This version is being made available in accordance with publisher policies. See http://orca.cf.ac.uk/policies.html for usage policies. Copyright and moral rights for publications made available in ORCA are retained by the copyright holders.



## An encapsulated Droplet Interface Bilayer Array for the High-Throughput Optical Measurement of Lipid Membranes with Single Bilayer Resolution.

Divesh K. Baxani, William D. Jamieson, David A. Barrow, Oliver K. Castell Cardiff University, Cardiff, United Kingdom

Droplet interface bilayers (DIBs) represent a recently developed method of generating stable artificial lipid membranes suitable for electrophysiology and optical interrogation. Prepared from the contacting of two aqueous droplets in oil, in the presence of dissolved lipid, DIBs are usually confined to a hydrophobic medium. Recent advances have created freestanding, hydrogel-encapsulated droplet interface bilayers (eDIBs) that are aqueous compatible, self-supporting and can withstand mechanical handling. Microfluidic methods allow for their rapid generation, paving the way for high throughput measurements on individually addressable lipid bilayers.

Here we report on the mass preparation of eDIBs using a 3D-printed microfluidic device and the proof-of-principle demonstration of high-throughput optical membrane screening with single bilayer resolution. Individual eDIBs are produced and output into wells of a 96-well plate. Optical measurements are made reporting on membrane leakage to monitor membrane integrity by fluorescent measurement. Addition of bilayer disrupting agents, such as detergents, to individual wells enables parallel measurement of membrane activity using a standard fluorescent plate reader This platform provides sufficient sensitivity to measure the leakage of dye through membrane spanning pores in otherwise intact bilayers, as well as direct detection of bilayer failure.

This high-throughput, scalable and automatable approach to arrayed bilayer measurements offers exciting opportunities for application with further fluorescentbased membrane, or protein assays. The segregated nature of each droplet allows for the separation of contents either side of the bilayer. Individual bilayers can be independently addressed, affording the opportunity to screen a wide parameter space with single bilayer resolution. The ability to reconstitute membrane proteins into DIBs and eDIBs creates many exciting opportunities for the optical screening of membrane proteins and high throughput biophysical measurements.